# Instruction for Use

- 1 **[Family name]** Ligature Directors  
 2 **[Manufacturer]** YDM CORPORATION  
 3 **[Model No and Device trade name]**

Model No.	Device trade name	UDI-DI	Material
22704	Ligature Director YS-704 Single End	4963931227044	(Working end) Stainless steel
22754	Ligature Director YS-704A Twister	4963931227549	
22755	Ligature Director YS-704B Double End	4963931227556	
27707	Ligature Director YS-704B AO	4963931277070	
27709	Arch Wire Director	4963931277094	(Handle) Brass
27721	Ligature Director YS-704B AO No Marking	4963931277216	
28002	Ligature Director Single End 3M Unitek	4963931280025	
28003	Ligature Director Double End 3M Unitek	4963931280032	

4 **[The device's intended purpose]**

Used to tuck and direct ligatures under the archwire or bracket wings or to push archwires or auxiliaries into position, also to ligate wires at "Ligature Director YS-704A Twister".

5 **[The performance characteristics of the device]**

<Overview>

This device has a variety of shapes to tuck and direct ligatures.

<Material>

Refer to "3 [Model No and Device trade name]"

6 **[How to use]**

Model No.	Device trade name	How to use
22704	Ligature Director YS-704 Single End	Tuck and direct ligatures under the archwire or bracket wings or push archwires or auxiliaries into position
22754	Ligature Director YS-704A Twister	One side, twists the device and ligature the wires. The other side, tucks and directs ligatures under the archwire or bracket wings or pushes archwires or auxiliaries into position.
22755	Ligature Director YS-704B Double End	Tuck and direct ligatures under the archwire or bracket wings or push archwires or auxiliaries into position
27707	Ligature Director YS-704B AO	
27709	Arch Wire Director	
27721	Ligature Director YS-704B AO No Marking	
22749	Cheek Retractor Doublewider S (2 ea.)	
22750	Cheek Retractor Doublewider L (2 ea.)	

7 **[General precautions and prohibition]**

- 1) The use of the product is strictly limited to orthodontic treatments. Precautions must be followed to ensure the appropriate use of the product.
- 2) Do not use the product for anything other than its intended use. The use of the product is limited to dental professionals. The use of the product is limited to dental professionals.
- 3) Properly clean, lubricate and sterilize before sterilization using the method and conditions specified by IFU for each patient.
- 4) Patients who have developed hypersensitivity symptoms such as rash and dermatitis or allergic symptoms due to treatment with this device should be discontinued and have a doctor diagnosed.
- 5) Do not change the design or remodel the product.
- 6) Do not drop the product. Avoid any strong impact on the product.
- 7) When disposing of the product, follow the instructions detailed in appropriate national regulations.
- 8) Do not heat product.
- 9) If chemicals adhere to the product, wipe it off immediately as it may corrode.

- 10) Stainless steel is more resistant to rust than iron, however it may corrode due to method of use and environment like temperature, humidity.

**8 [Limitations on processing]**

Stop using the product immediately if deterioration or abnormality are visible.

**9 [Initial treatment at the point of use]**

Properly clean, lubricate and sterilize before sterilization using the method and conditions specified by IFU for each patient.

**# [Preparation before cleaning]**

See next section.

**11 [The process of cleaning and sterilization/disinfection]**

**1) General cautions**

- a) When processing this product, please refer to the "The Fundamental Guide Book for Cleaning and Sterilizing Instruments" published by our company or the item of "maintenance" on our website (<http://www.ydm.co.jp/>).
- b) The device contaminated with blood, body fluid, tissue fragments, chemicals, etc. shall be cleaned and sterilized immediately in compliance with following method in order to prevent The contaminants from drying out and sticking. Leaving The contaminants attached can make it difficult to remove.

**2) Procedure**

- a) Clean the device with a cleaning agent in the method of i ) or ii ).
  - i ) Using an ultrasonic cleaner with medical rust preventive cleaner [for example, "Z-1 eco (sold separately)"] is dissolved in warm water, and ultrasonic cleaning is performed for 5 minutes.
  - ii ) Immerse and wash in a rust preventive detergent for medical use for 15 minutes.
- b) Wash off dirt and cleaning liquid adhered to the device with running water. (If soap components remain in the cleaning solution and remain attached, it may cause rust and stains.)
- c) Allow the device to dry. (Remaining moisture may cause rust and reduced sterilization effect)
- d) Sterilize using an autoclave sterilizer.

Conditions: 134°C to 138°C, 3 minutes to 3.5 minutes

\* To prevent scratches due to contact with other devices, it is recommended to put them in a sterilized pack or case.

e) After use, Clean immediately according to a) above.

**3) Cautions of cleaning and sterilization/disinfection**

- a) Do not use the following chemicals as they may cause metal corrosion.  
sodium hypochlorite, benzalkonium chloride, chlorhexidine gluconate, tincture of iodine, Iodoform and peracetic acid may corrode metal.
- b) Do not perform unsaturated Chemical Vapor (Chemical Autoclave) on the product, as such agents may cause corrosion.
- c) Do not perform plasma sterilization because it affects the material.
- d) Do not use household cleaning agents as such chemical may lead to corrosion.  
Use of Antirust Cleaner(13-538 Z-1 eco Fine Liquid) for dental instrument is recommended.
- e) Wire brush and steel wool may cause some damage for instruments.  
Do not use such materials when you clean instruments
- f) Do not use Super Acidic Water as such agents may cause corrosion.  
Purified water (no impurities ) is recommended for rinsing.  
Tap water contains chlorine. This chlorine may corrode the metal.
- g) When using a cleaning device such as a washer-disinfector, the attached dirt or cleaning solution may not be completely removed if water saving or a time saving program is used. Be sure to refer to the manufacturer's instruction manual, etc. to ensure rinsing and to remove dirt.

- h) Note the following when handling autoclave sterilizers.
- i ) The autoclave sterilizer may be contaminated in the cabinet depending on the usage status and period. If autoclave sterilization is performed with dirt attached, stains may adhere to the instrument. Encourage regular cleaning according to the package insert or instruction manual of the sterilizer so that the inside of the chamber does not become dirty. In particular, the chamber lid packing and air filter may require periodic replacement.
  - ii ) Use purified water (pure water) as much as possible. When tap water is used, the equipment may corrode due to the influence of chlorine ions.
  - iii ) Pay attention to the drying temperature and internal temperature.
  - iv ) Do not place this product near the heater. (It may be higher than the indicated temperature)
  - v ) If there is a risk of the inside of the chamber becoming hot, preheat drying. High temperature drying may cause the instrument to be altered or discolored, deteriorated or damaged.
  - vi ) If autoclave sterilization is performed while cleaning and rinsing are not complete, or when water is adhering to the inside of the autoclave sterilizer chamber, there is a risk of seizure of the instrument.
- i) After washing and sterilization, remove the water adhering to this instrument and dry it thoroughly before storing. Leaving it for a long time with moisture attached may cause rusting, spots, etc.

## 12 [Maintenance, Inspection and Testing]

After cleaning and sterilization, check the following before use.

- 1) Check that there is no dirt, damage, cracks, scratches.

## 13 [Storage]

- 1) Store in a clean place free from dust and dirt, avoiding high temperature and humidity.
- 2) To prevent galvanic corrosion due to metal potential difference, do not store instruments of different materials together.
- 3) Note the following in order to prevent rust stains.
  - a) Do not store with rusted instruments.
  - b) Do not store with chemicals.
  - c) Pay attention to rust generated inside sterilizers and storages.
- 4) Carefully move and store in a tray / container as this may cause deformation. When using trays or containers, do not place heavy equipment.

## 14 [Disposal considerations]

Dental specialists and / or orthodontists are responsible for disposal in accordance with relevant applicable law.

YDM CORPORATION by improper disposal. expressly disclaims any liability for the spread of illness or

## 15 [Additional information]

<Breakage or Sharpening>

In case of breakage of product, please report to dealer from which you purchased in format prepared by YDM CORPORATION named "Claim Report".

Please see more details subject "K".

In case of sharpening of cutters, we will do with a certain of charges.

## 16 [Manufacturer contact]

**YDM CORPORATION**

28 Imaizumi, Higashi-Matsuyama-shi,

Saitama 355-0042, Japan

TEL: 81-3-3828-3161 FAX: 81-3-3827-8991

URL: <http://www.ydm.co.jp/> E-mail: [ydm@ydm.co.jp](mailto:ydm@ydm.co.jp)

17 **[Notice to the user and/or patient]**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

18 **[Additional information]**

In case of defect of device, please report to dealer from which you purchased in format prepared by YDM Corporation named "Claim Report".

19 **[Description of symbols]**

Medical device

Manufacturere

Lot number

Notice

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Ver.01 7th September, 2021

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1 **[Familiennamen]** Ligature Directors

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<Materyal>

Bkz. "3 [Model No ve Cihazın Ticari Adı]"

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Model No.	Cihazın ticari adı	Kullanımı
22704	Ligature Director YS-704 Single End	Tuck and direct ligatures under the archwire or bracket wings or push archwires or auxiliaries into position
22754	Ligature Director YS-704A Twister	One side, twists the device and ligature the wires. The other side, tucks and directs ligatures under the archwire or bracket wings or pushes archwires or auxiliaries into position.
22755	Ligature Director YS-704B Double End	Tuck and direct ligatures under the archwire or bracket wings or push archwires or auxiliaries into position
27707	Ligature Director YS-704B AO	
27709	Arch Wire Director	
27721	Ligature Director YS-704B AO No Marking	
22749	Cheek Retractor Doublewider S (2 ea.)	
22750	Cheek Retractor Doublewider L (2 ea.)	

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