Instruction for Use

- 1 **[Family name]** Photographic Mirror
 2 **[Manufacturer]** YDM CORPORATION
 3 **[Model No and Device trade name]**

| Model No. | Device trade name | UDI-DI | Material | |
|-----------|--|---------------|-----------------|-----------------------|
| 13261 | Photographic Mirror ST for Occlusal | 4963931132614 | Stainless steel | |
| 13262 | Photographic Mirror ST for Lingual | 4963931132621 | | |
| 13263 | Photographic Mirror ST for Buccal | 4963931132638 | | |
| 13264 | Photographic Mirror ST for Occlusal Wide L | 4963931132645 | | |
| 13265 | Photographic Mirror ST for Occlusal Wide S | 4963931132652 | | |
| 13266 | Photographic Mirror ST for Occlusal Wide LL | 4963931132669 | | |
| 13271 | Photographic Mirror SE for Buccal & Lingual | 4963931132713 | | |
| 13272 | Photographic Mirror SE for Occlusal L | 4963931132720 | | |
| 13273 | Photographic Mirror SE for Occlusal M | 4963931132737 | | |
| 13274 | Photographic Mirror SE for Occlusal LL | 4963931132744 | | |
| 13281 | Photographic Mirror Crystal for Buccal & Lingual | 4963931132812 | | Glass (Metal film) |
| 13282 | Photographic Mirror Crystal for Occlusal L | 4963931132829 | | |
| 13283 | Photographic Mirror Crystal for Occlusal M | 4963931132836 | | |
| 13571 | Photographic Mirror for Occlusal | 4963931135714 | Stainless steel | |
| 13572 | Photographic Mirror for Lingual | 4963931135721 | | |
| 13573 | Photographic Mirror for Buccal | 4963931135738 | | |
| 13576 | Photographic Mirror for Occlusal Wide L | 4963931135769 | | |
| 13577 | Photographic Mirror for Occlusal Wide S | 4963931135776 | | |
| 13578 | Photographic Mirror for Occlusal Wide LL | 4963931135783 | | |
| 13641 | Contrafter for Anterior L | 4963931136414 | Aluminum | |
| 13642 | Contrafter for Posterior | 4963931136421 | | |
| 13643 | Contrafter for Anterior S | 4963931136438 | | |

4 **[The device's intended purpose]**
 Used to assist practitioners during general surgery by displaying the image of an object placed before it.

5 **[The performance characteristics of the device]**
 <Overview>
 This device has a variety of shapes to assist practitioners during general surgery by displaying the image of an object placed before it.
 <Material>
 Refer to "3 [Model No and Device trade name]"

6 **[How to use]**
 Insert it into the oral cavity and display the state of the oral cavity.

- 7 **[General precautions and prohibition]**
- 1) The use of the product is strictly limited to dental treatments.
 Precautions must be followed to ensure the appropriate use of the product.
 - 2) Do not use the product for anything other than its intended use. The use of the product is limited to dental professionals.
 The use of the product is limited to dental professionals.
 - 3) Properly clean, lubricate and sterilize before sterilization using the method and conditions specified by IFU for each patient.
 - 4) Patients who have developed hypersensitivity symptoms such as rash and dermatitis or allergic symptoms due to treatment with this device should be discontinued and have a doctor diagnosed.
 - 5) Do not change the design or remodel the product.
 - 6) Do not drop the product. Avoid any strong impact on the product.

- 7) When disposing of the product, follow the instructions detailed in appropriate national regulations.
- 8) Do not heat product.
- 9) If chemicals adhere to the product, wipe it off immediately as it may corrode.
- 10) Stainless steel is more resistant to rust than iron, however it may corrode due to method of use and environment like temperature, humidity.

8 [Limitations on processing]

Stop using the product immediately if deterioration or abnormality are visible.

9 [Initial treatment at the point of use]

Properly clean, lubricate and sterilize before sterilization using the method and conditions specified by IFU for each patient.

10 [Preparation before cleaning]

See next section.

11 [The process of cleaning and sterilization/disinfection]

1) General cautions

- a) When crossing this product, please refer to the "The Fundamental Guide Book for Cleaning and Sterilizing Instruments" published by our company or the item of "maintenance" on our website (<http://www.ydm.co.jp/>).
- b) The device contaminated with blood, body fluid, tissue fragments, chemicals, etc. shall be cleaned and sterilized immediately in compliance with following method in order to prevent The contaminants from drying out and sticking. Leaving The contaminants attached can make it difficult to remove.

2) Procedure

"Photographic Mirror"

- a) Immerse and clean in a rust preventive detergent for medical use for 15 minutes.
Shall not use an ultrasonic cleaner as it may damage the mirror surface.
- b) Wash off dirt and cleaning liquid adhered to the device with running water (If soap components remain in the cleaning solution and remain attached, it may cause rust and stains.)
- c) Use a soft cloth such as a microfiber cloth that does not generate lint, and remove the water without rubbing.
- d) Put in the attached protective bag.
- e) Enclose in a sterilization bag. (If it is sealed directly in a sterilization bag without using a protective bag, stains due to the ink and adhesive components of the sterilization bag may occur.)
- f) Place the Photographic Mirror into the autoclave sterilizer with the mirror side facing down and the transparent film side of the sterilization bag facing up.
(This allows to prevent stains caused by water droplets generated during vacuum cooling of autoclave sterilization.)
- g) Sterilize using an autoclave sterilizer.
Conditions: 134°C to 138°C, 3 minutes to 3.5 minutes
- h) After use, Clean immediately according to a) above.

"Contraster"

- a) Clean the device with a cleaning agent in the method of i) or ii).
 - i) Using an ultrasonic cleaner with medical rust preventive cleaner [for example, "Z-1 eco (sold separately)"] is dissolved in warm water, and ultrasonic cleaning is performed for 5 minutes.
 - ii) Immerse and wash in a rust preventive detergent for medical use for 15 minutes.
 - b) Wash off dirt and cleaning liquid adhered to the device with running water. (If soap components remain in the cleaning solution and remain attached, it may cause rust and stains.)
 - c) Allow the device to dry. (Remaining moisture may cause rust and reduced sterilization effect)
 - d) Sterilize using an autoclave sterilizer.
Conditions: 134°C to 138°C, 3 minutes to 3.5 minutes
- * To prevent scratches due to contact with other devices, it is recommended to put them in a sterilized pack or case.
- e) After use, Clean immediately according to a) above.

- 3) Cautions of cleaning and sterilization/disinfection
- a) Do not use the following chemicals as they may cause metal corrosion.
sodium hypochlorite, benzalkonium chloride, chlorhexidine gluconate, tincture of iodine, Iodoform and peracetic acid may corrode metal.
 - b) Do not perform unsaturated Chemical Vapor (Chemical Autoclave) on the product, as such agents may cause corrosion.
 - c) Do not perform plasma sterilization because it affects the material.
 - d) Do not use household cleaning agents as such chemical may lead to corrosion.
Use of Antitrust Cleaner(13-538 Z-1 eco Fine Liquid) for dental instrument is recommended.
 - e) Wire brush and steel wool may cause some damage for instruments.
Do not use such materials when you clean instruments
 - f) Do not use Super Acidic Water as such agents may cause corrosion.
Purified water (no impurities) is recommended for rinsing.
Tap water contains chlorine. This chlorine may corrode the metal.
 - g) When using a cleaning device such as a washer-disinfector, the attached dirt or cleaning solution may not be completely removed if water saving or a time saving program is used. Be sure to refer to the manufacturer's instruction manual, etc. to ensure rinsing and to remove dirt.
 - h) Note the following when handling autoclave sterilizers.
 - i) The autoclave sterilizer may be contaminated in the cabinet depending on the usage status and period. If autoclave sterilization is performed with dirt attached, stains may adhere to the instrument. Encourage regular cleaning according to the package insert or instruction manual of the sterilizer so that the inside of the chamber does not become dirty. In particular, the chamber lid packing and air filter may require periodic replacement.
 - ii) Use purified water (pure water) as much as possible. When tap water is used, the equipment may corrode due to the influence of chlorine ions.
 - iii) Pay attention to the drying temperature and internal temperature.
 - iv) Do not place this product near the heater. (It may be higher than the indicated temperature)
 - v) If there is a risk of the inside of the chamber becoming hot, preheat drying. High temperature drying may cause the instrument to be altered or discolored, deteriorated or damaged.
 - vi) If autoclave sterilization is performed while cleaning and rinsing are not complete, or when water is adhering to the inside of the autoclave sterilizer chamber, there is a risk of seizure of the instrument.
 - i) After washing and sterilization, remove the water adhering to this instrument and dry it thoroughly before storing. Leaving it for a long time with moisture attached may cause rusting, spots, etc.

12 [Maintenance, Inspection and Testing]

After cleaning and sterilization, check the following before use.

- 1) Check that there is no dirt, damage, cracks, scratches.

13 [Storage]

- 1) Store in a clean place free from dust and dust, avoiding high temperature and humidity.
- 2) To prevent galvanic corrosion due to metal potential difference, do not store instruments of different materials together.
- 3) Note the following in order to prevent rust stains.
 - a) Do not store with rusted instruments.
 - b) Do not store with chemicals.
 - c) Pay attention to rust generated inside sterilizers and storages.
- 4) Carefully move and store in a tray / container as this may cause deformation. When using trays or containers, do not place heavy equipment.

14 **[Disposal considerations]**

Dental specialists and / or orthodontists are responsible for disposal in accordance with relevant applicable law. YDM CORPORATION by improper disposal. expressly disclaims any liability for the spread of illness or personal injury caused

15 **[Additional information]**

<Breakage or Sharpening>

In case of breakage of product, please report to dealer from which you purchased in format prepared by YDM CORPORATION named "Claim Report".

Please see more details subject "K".

In case of sharpening of cutters, we will do with a certain of charges.

16 **[Manufacturer contact]**

YDM CORPORATION

28 Imaizumi, Higashi-Matsuyama-shi,

Saitama 355-0042, Japan

TEL: 81-3-3828-3161 FAX: 81-3-3827-8991

URL: <http://www.ydm.co.jp/> E-mail: ydm@ydm.co.jp

17 **[Notice to the user and/or patient]**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

18 **[Additional information]**

In case of defect of device, please report to dealer from which you purchased in format prepared by YDM Corporation named "Claim Report".

19 **[Description of symbols]**

Medical device

Manufacturere

Lot number

Notice

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

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- 1 **[Familiename]** Photographic Mirror
 2 **[Hersteller]** YDM CORPORATION
 3 **[Modellnr. und Handelsname]**

| Modellnr. | Handelsname | UDI-DI | Material |
|-----------|--|---------------|-------------------------|
| 13261 | Photographic Mirror ST for Occlusal | 4963931132614 | Edelstahl |
| 13262 | Photographic Mirror ST for Lingual | 4963931132621 | |
| 13263 | Photographic Mirror ST for Buccal | 4963931132638 | |
| 13264 | Photographic Mirror ST for Occlusal Wide L | 4963931132645 | |
| 13265 | Photographic Mirror ST for Occlusal Wide S | 4963931132652 | |
| 13266 | Photographic Mirror ST for Occlusal Wide LL | 4963931132669 | |
| 13271 | Photographic Mirror SE for Buccal & Lingual | 4963931132713 | |
| 13272 | Photographic Mirror SE for Occlusal L | 4963931132720 | |
| 13273 | Photographic Mirror SE for Occlusal M | 4963931132737 | |
| 13274 | Photographic Mirror SE for Occlusal LL | 4963931132744 | |
| 13281 | Photographic Mirror Crystal for Buccal & Lingual | 4963931132812 | Glas (Metallschicht) |
| 13282 | Photographic Mirror Crystal for Occlusal L | 4963931132829 | |
| 13283 | Photographic Mirror Crystal for Occlusal M | 4963931132836 | |
| 13571 | Photographic Mirror for Occlusal | 4963931135714 | Edelstahl |
| 13572 | Photographic Mirror for Lingual | 4963931135721 | |
| 13573 | Photographic Mirror for Buccal | 4963931135738 | |
| 13576 | Photographic Mirror for Occlusal Wide L | 4963931135769 | |
| 13577 | Photographic Mirror for Occlusal Wide S | 4963931135776 | |
| 13578 | Photographic Mirror for Occlusal Wide LL | 4963931135783 | |
| 13641 | Contrafter for Anterior L | 4963931136414 | Aluminium |
| 13642 | Contrafter for Posterior | 4963931136421 | |
| 13643 | Contrafter for Anterior S | 4963931136438 | |

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